


ENDOJACK UTERINE MANIPULATOR

A sterile disposable kit for the movement of the uterus



The uterine manipulator is designed for the mobilisation of the uterus in the pelvic floor during laparoscopic surgical procedures



FEATURES

- A rotating handle to enable the movement of the tip inside the uterus
- An inflatable balloon for securing the device inside the uterus
- An adjustable spacer with two positions to define the penetration depth
- Two lines with luer fittings to inflate the balloon and to inject the contrast liquid



ADVANTAGES

- Ready to use device
- Ideal lightness
- Easy-to-use kit